easyScan™

Distribué par / Distributed by
www.visiomed-lab.com

Marché Européen / European Market
Visiomed Group SA
Service clients / Customer service
8, avenue Kléber
75116 Paris - France
Tel : +33 8 92 350 334
Fax : +33 1 43 73 59 56
contact@visiomed-lab.fr

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By **Visiomed**®

Dear Customer,

You have just purchased this **VISIONED® EasyScan™ VM-ZX1** temporal and auricular thermometer and we thank you. We wish you excellent use of the product and recommend you read this manual.

CONTENTS

1- SAFETY INSTRUCTIONS	20
2- PRECAUTIONS BEFORE USE	22
3- INTRODUCTION	23
4- BODY TEMPERATURE	23
5- MAIN FUNCTIONS	24
6- DESCRIPTION OF THE UNIT	25
7- LCD SCREEN DISPLAY	25
8- SUGGESTIONS BEFORE MEASURING	25
9- CHOICE OF TEMPERATURE UNIT	27
10- HOW TO MEASURE TEMPERATURE	27
11- DATA STORAGE	28
12- BATTERY REPLACEMENT	28
13- UPKEEP, STORAGE, AND CALIBRATION	29
14- BREAKDOWNS AND REPAIRS	29
15- SPECIFICATIONS	30

THE MANUFACTURER RESERVES THE RIGHT TO MODIFY THE
PRODUCT SPECIFICATIONS WITHOUT PRIOR NOTICE

1- SAFETY INSTRUCTIONS

- Follow the maintenance instructions specified in this manual.
- Use this unit only for its intended purpose as described in this manual.
- This unit can be used for personal use at home.
- Use this unit within an ambient temperature range between 15° and 39.9°C.

- This unit should always be stored in a clean, dry place.
 - Do not expose this unit to extreme temperature conditions $T^{\circ} > 55^{\circ}\text{C}$ or $T^{\circ} < -20^{\circ}\text{C}$.
 - Do not use this unit at a relative humidity above 85%.
 - Do not expose this unit to sunlight or water.
 - Do not use this unit outdoors.
 - Do not expose this unit to electrical shocks.
 - Never drop this unit.
 - Do not perform maintenance on this unit yourself. Do not attempt to open the unit. In the event of a problem, contact your retailer.
 - Do not leave this unit within the reach of children. A self-measurement taken forcibly by children can damage their ear. In the event of accidental ingestion of the battery or the protective film, see a physician immediately.
 - The measurement results are for informational purposes only. If unsure about your results, please contact your physician.
 - This product should not encourage self-medication or changes to your treatment.
 - Stop using the unit in the event of a failure or malfunction.
 - This electrical medical equipment requires special precautions regarding electromagnetic compatibility. It must be installed and used in accordance with electromagnetic information provided in the EMC Data table.
 - Readings may be distorted if the unit is used near a television, microwave oven, mobile phone or any other device with an electrical field.
 - Do not use battery other than the battery mentioned, do not recharge non-rechargeable battery, do not incinerate them.
 - Remove the battery when not using the thermometer over a long period.
 - Respect the polarity of the battery. A polarity error can cause damage and void the warranty of your unit.
 - Stop using the unit if any pain occurs*. There is a risk of damaging the ear canal.
 - It is not recommended to use the device on patients with ear diseases such as otitis externa and tympanitis.
- *The affected area may suffer damage.*
- Do not use the unit in a wet ear canal, such as after swimming or after a

bath. There is a risk of damaging the ear canal.

- Do not touch and or blow on the infrared sensor. A dirty or damaged infrared sensor may cause inaccurate results.

WARNING: You cannot dispose of this unit with your other household waste. You are required to dispose of the unit at a designated collection point for recycling e-waste. For more information on e-waste collection points, please contact your town hall or your household waste disposal service.

2- PRECAUTIONS BEFORE USE

The **Easyscan™ VM-ZX1** thermometer is pre-configured at the factory. It is not necessary to calibrate the device when it is first used.

In order to obtain reliable, stable results, it is advisable, with each change of environment with a significant difference in ambient temperature, to leave the **Easyscan™ VM-ZX1** at this ambient temperature for 15 to 20 minutes before use. It is important to wait a period of 1 minute between readings.

For a temporal temperature reading (right temporal artery) take care to brush away the hair, wipe any sweat off the brow, avoid air flows, (e.g. supplemental oxygen, air conditioning, etc.).

To measure the tympanic temperature, make sure that the ear canal is clean, as an ear with excess wax can give inaccurate readings.

The unit cannot be reused if the screen is not turned off.

- Do not place the thermometer in contact with diseased, ulcerated, or bruised skin.
- Do not drink hot or cold beverages, and do not perform strenuous exercise during the temperature reading.
- Watch out for physiological temperature changes that are to be accounted for depending on the results: Body temperature increases by 0.5°C between 6 am and 3 pm. Women have a body temperature about 0.2°C higher on average. Their body temperature also varies depending on the ovulation cycle. Thus it increases by 0.5°C in the second half of the cycle, and in early pregnancy.
- When seated, the body temperature is 0.3 to 0.4°C lower relative to the standing position.

3- INTRODUCTION

Easyscan™ VM-ZX1, developed by Visiomed®, is an Infrared electronic medical thermometer with contact equipped with the latest in automatic calibration technology: **MicroSecond Flash™**.

Easy to use, it enables temporal (right temporal artery) or ear temperature measurement, in less than a second. Fast and accurate, just put the **Easyscan™ VM-ZX1** probe against the right temple or in the ear for an instant temperature reading. **Easyscan™** is particularly effective for taking the temperature of infants and young children. It has shown that temperature readings taken on the temporal artery are, in new-borns, more accurate than tympanic thermometry and better tolerated than rectal thermometry*. **Easyscan™** is suitable for the whole family.

* *Greenes D, Fleisher G. Accuracy of a Non-invasive Temporal Artery Thermometer for Use in Infants. Arch Pediatr Adolesc Med 2001; 155:376.*

Its control LEDs, green, orange or red, indicate the temperature result in colour, for a better diagnosis.

This product should not be used for any purpose other than temporal (right temporal artery) or ear temperature measurement. It can be used on all persons, new-borns, children or adults.

4- BODY TEMPERATURE

4.1. IMPORTANT CONCEPT:

ESTIMATING BODY TEMPERATURE

The normal body temperature is located within a range, and a fixed value. Everyone has their own temperature range.

Body temperature can vary by age, gender, time of day, activity, and emotions. Persons of different ages have different temperatures; likewise, body temperature varies at different times of the day. For example, children's temperatures are 0.5°C higher than those of adults; most people have a lower body temperature at night, and it increases throughout the day. Different parts of the body yield give different results. Normally, the rectal temperature is 0.3°C higher than the oral temperature, and the oral temperature is 0.3°C higher than the temperature of the armpit.

For your own temperature and that of your family members, you can measure the temperature at different times of the day over 3 or 4 days, when you are healthy. Once you know the temperature range of each

person, you can save it and compare it to the temperature when you are measuring for fever.

4.2. BODY TEMPERATURE REGULATION

When the body temperature rises or falls, the body initially regulates its temperature from the brain. Since the temple and the ear are the body parts nearest the brain, they are the first to be informed of a temperature rise.

4.3. HOW TO PROPERLY USE THE THERMOMETER IN TEMPORAL MODE

SEE DIAGRAM **A** PAGE 3.

Easyscan™ VM-ZX1 is pre-configured at the factory. It is not necessary to calibrate the device when it is first used.

Keep the tip containing the temporal probe on the unit.

Place the thermometer against the right temple, where the right temporal artery is located.

4.4. HOW TO PROPERLY USE THE THERMOMETER IN TYMPANIC (EAR) MODE

SEE DIAGRAM **A** PAGE 3.

Remove the temporal probe tip from the unit. For a child under one year of age: pull the ear diagonally to the rear, to straighten the ear canal, and gently insert the probe to take the reading. For a child over one year of age, and adults: pull the ear back slightly and gently insert the probe along the ear canal to take the reading.

5- MAIN FUNCTIONS

1. Display of the temperature reading in degrees Celsius or Fahrenheit: Temperature range from 32°C to 42.9°C (89.6°F to 109.2°F).
2. Memorization of the last 30 readings taken: The temperature is recorded.
3. Audible alert during temperature measurement based on the result.
4. Back-lit colour display to indicate and quickly interpret the diagnosis of the temperature reading.

6- DESCRIPTION OF THE UNIT

SEE DIAGRAM **B** PAGE 3.

- | | |
|------------------------------|------------------------------|
| ① Temporal measurement probe | ② Tympanic measurement probe |
| ③ LCD Screen | ④ Measurement |
| ⑤ °C/°F button | ⑥ MEM Button |
| ⑦ Battery | |

7- LCD SCREEN DISPLAY

SEE DIAGRAM **C** PAGE 3.

- ① Tympanic Mode / Temporal Mode
- ② Ambient air
- ③ Viewing the measurement mode.

8- SUGGESTIONS BEFORE MEASURING

1. Replace batteries that are too drained to effectively power the unit.
2. For a proper reading, condition the thermometer by placing it in the environment where the measurement will be taken 30 minutes before use
3. An unexpected fluctuation of the ambient temperature may distort the accuracy of the reading. For example, if the temperature of the place where the measurement is taken is different from the temperature where the thermometer is stored, when one wants to take a temperature reading near an air conditioner. The results will be unreliable.
4. The measurement should be taken at rest, to yield accurate results. Body temperature is likely to increase just after exercise or after a bath.
5. Before taking a temperature, in order to obtain a stable, reliable result
 - 5.1 Temporal mode:
 - Take care to brush away the hair.
 - Wipe any sweat from the brow.
 - Avoid air flows (e.g. supplemental oxygen, air conditioning, etc.).
 - Wait a period of 1 minute between readings.
 - With each change of environment having a high ambient temperature difference, leave the thermometer at this ambient temperature for at least 30 min before use.

5.2 Tympanic mode:

The temperature in the right ear may differ from that of the left ear. Therefore, you should always take the temperature in both ears and use the highest temperature as a reference. The ear must be free of excessive wax build-up or plugs, so that the temperature measurement is accurate. In the event of otitis, **Easyscan™** should be used in temporal mode. External factors can influence the temperature reading in the ear, as in the following cases where the patient:

- Is lying on one ear
- Has their ears covered
- Is exposed to very hot or very cold temperatures
- Recently went swimming or took a bath.

In the above cases, place the patient in another environment and wait 20 minutes before taking their temperature. If the event of auricular therapy (prescription ear drops or other internal medication), be careful not to take the ear temperature of the treated ear. Since the ear canal is under treatment, the temperature reading displayed by **Easyscan™** will not reflect the actual body temperature.

Recommendations

1. For inexperienced users, it is recommended to measure the same ear 3 times during the first use. The maximum value will be saved if different results are recorded, for every thermometer under continuous operation may show a slight margin of error (+/- 0.3°C, +/- 0.5°F).
2. Wait about 60 seconds before measuring the same patient again to avoid excessive cooling of the skin.
3. For healthy users, both ears are equivalent in terms of body temperature. Choose the ear that has a relatively stable temperature, higher than the other ear.
4. Please clean the ear before measurement

9- CHOICE OF TEMPERATURE UNIT

Press the °C/°F button to select the temperature unit. The display shows °C to display degrees Celsius or °F to display degrees Fahrenheit.

10- HOW TO MEASURE TEMPERATURE

Insert the batteries. For the first use or when inserting new batteries, wait about 10 minutes to preheat the device.

10.1 For a temporal measurement: Keep the tip containing the temporal probe on the unit. Place the probe against the right temple (SEE DIAGRAM A PAGE 3). Then press the button once to start the temperature reading and release it: the unit will beep once. When the reading is complete the unit will beep 2 times. The result is displayed on the screen for one minute then the ambient temperature appears on the screen for one minute too, with the logo .

10.2 For a tympanic temperature measurement: Remove the tip of the temporal probe from the unit, insert the tympanic probe along the ear canal. (SEE DIAGRAM A PAGE 3) Then press the button once to start the temperature reading and release it: the unit will beep once. When the reading is complete, the unit will beep 2 times. The result is displayed on the screen.

- Visual diagnosis: Back-lit colour display to indicate and quickly interpret the diagnosis of the temperature reading.
- Audio diagnosis:
 - When the reading is less than 37.5°C, the unit will make 3 short beeps.
 - When the reading is higher than 37.5°C, the unit will make 3 triple beeps.

Before any temperature measurements, take care to brush away the hair and clear any sweat from the forehead.

The result is automatically stored. When taking consecutive temperature measurements, it is possible that the results will not be identical. Wait at least 1 minute between each reading.

Displaying temperature results:

1. Measurement range: 32°C - 42.9°C (89.6°F - 109.2°F)
2. *HI* is displayed for a result above 42.9°C (109.2°F).
3. *LO* is displayed for a result below 32°C (89.6°F).

Automatic shutdown of the back-lighting after 8 ~ 15 seconds.

11- DATA STORAGE

Open the battery compartment to access the settings buttons.

Press the *MEM* button, the LCD screen reads  and displays the last temperature reading.

1. *--* will be displayed if no data is stored
2. Hit the *MEM* button to scroll through the previous data, then hit the *°C/°F* to scroll through the next data.

Note: When changing batteries, the stored data will be erased.

12- BATTERY REPLACEMENT

When the screen reads , the batteries are drained.

Remove the battery compartment cover and replace the old battery, taking care to respect the polarity.

The LCD screen will display all the symbols on the screen.

Use only new batteries (1 x CR2032 button battery), never use rechargeable batteries. Do not mix an old battery with a new battery.

Replace the cover after inserting the battery.

IMPORTANT : Respect the polarity of the batteries. A polarity error can cause damage and void the warranty of your unit.

13- UPKEEP, STORAGE, AND CALIBRATION

UPKEEP:

1. Remove any spots on the case with a soft, dry cloth.
2. Clean the probe tip as follows: wipe the surface very gently with a buffer, a soft cloth lightly moistened with 70% surgical spirit; you can use the thermometer once the alcohol has completely evaporated.
3. Do not immerse the unit in any liquid whatsoever.
4. When dirty, gently clean the IR sensor with a soft, dry cloth.

Do not clean the infrared sensor with tissue or paper towel that may scratch it, which could lead to inaccurate results.

STORAGE:

1. Do not expose the thermometer to sunlight, high temperatures, or under conditions of high humidity, dust, near fire; do not subject it to vibrations or direct impacts.
2. Remove the batteries from the housing if you do not plan to use the unit over a long period.

CALIBRATION:

This thermometer is pre-configured at the factory. If the operating instructions are followed, a re-calibration is not necessary. If you have doubts about the accuracy of temperature readings, please contact your retailer.

Do not disassemble the thermometer yourself - this will void your warranty.

14- BREAKDOWNS AND REPAIRS

- **No response:** Check the battery (polarity, wear)
- **The battery symbol appears:** Change the batteries. Check battery polarity during insertion. A polarity error can damage your unit.
- **The symbol  symbol appears:** Wrong measurement
- **Lo appears on screen** (result below 32°C or 89.6°F): Is the thermometer placed properly against the temple or in the ear? Follow the instructions in the user manual regarding the application of the probe according to the different measurement modes.
- **Hi appears on screen** (result above 42.9°C or 109.2°F): Please check the procedure. Follow the instructions in the operating manual to take readings properly.

15- SPECIFICATIONS

Type of Product	Infrared medical thermometer
Model	EasyScan™ VM-ZX1
Medical device	Class IIa
Temperature sensor	IR Sensor
Wavelength	5-14µm
Emissivity	0.95
Normal conditions of use	Auricular : 10-39.9°C (50-103.9°F) Temporal : 15-39.9°C (59-103.9°F)
Storage conditions	Temperature: -20°C - 55°C (-4 - 131°F) Humidity: 30% RH ~ 85% RH
Atmospheric pressure	101.325Kpa
Humidity level	30% RH ~ 80% RH
Power	1 x CR2032 button battery
Dimensions	140,3 x 45,7 x 43,8mm (±2mm)
Weight	54,6g (without battery)
Displaying of readings	°C / °F
Body temperature range	32°C - 42.9°C (89.6°F - 109.2°F)
Repeatability	≤±0.3°C
Display resolution	0.1°C
Accuracy	35.5°C ~ 42°C : ±0.2°C <35.5°C : ±0.3°C, >42°C : ±0.3°C
Data storage	Last 30 readings taken
Consumption	<75 mW
Voltage	3V ~ 6V
Idle current	≤20µA
Dynamic current	≤25mA
Automatic shutoff	60 seconds

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EMC DATA

EMC data

IMPORTANT INFORMATION REGARDING ELECTROMAGNETIC COMPATIBILITY (EMC)

With the increased number of electronic devices such as computers and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation. Medical devices should also not interfere with other devices. In order to regulate the requirements for EMC (Electromagnetic Compatibility) with the aim to prevent unsafe product situations, the IEC60601- 1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices. The devices conform to this IEC60601-1-2:2007 standard for both immunity and emissions.

The device needs to be installed and put into service in accordance with the information provided in the user manual.

Nevertheless, special precautions need to be observed:

- The use of accessories other than those provided by manufacturer, may result in increased emission or decreased immunity of the device.
- Refer to EMC table guidance regarding the EMC environment in which the device should be used.

Guidance and manufacturer's declaration - electromagnetic immunity. The device is intended for use in the electromagnetic environment specified below. The customer or the user assures that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration - electromagnetic immunity. The device is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.

Immunity test	IEC 60601	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

Note : UT is the a.c mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity. The device is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.

Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	$d = 1,2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2.5 GHz
			where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). b Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V/m].

Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance *d* away from the equipment. The distance *d* is calculated by the MANUFACTURER from the 800 MHz to 2,5 GHz column of Table 6 of IEC 60601-1-2:2007, as appropriate.

Recommended separation distances between portable and mobile RF communications equipment and the device .

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile rf communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter	Separation distance according to frequency of transmitter m		
	150 KHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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EN: Visiomed® will repair or replace this product free of charge in the case of defective parts or manufacturing defects, in accordance with the conditions mentioned below as follows:

DURATION: 24 MONTHS RETURN TO WORKSHOP

LIMITS AND EXCLUSIONS: This guarantee concerns only the original final purchaser. A purchase invoice, or another proof of purchase, with this guarantee card will be required to obtain an after-sales service, in accordance with this guarantee. This guarantee card will not be extended to another person only the original final purchaser. This guarantee becomes void if the serial numbers on the product are modified, replaced, illegible, absent, or if repair has been carried out by a service not approved, including the user.

This guarantee covers only the defects of the material or parts, occurring during normal use of the product. It does not cover the damage caused during the transport of the apparatus, causes due to repairs being carried out by the distributor, by any modifications undertaken, any connection of equipment not approved by Visiomed®, or causes contrary to those written in the user manual or notice. Moreover, the present guarantee does not cover damage due to falls, bad handling, bad installations, damage by fire, floods, lightning, or any other natural disaster. This guarantee does not cover the packing of the material, the accessories, the defects caused by commercial exposure of the product, show room, sale space, demonstration etc.... Normal maintenance, cleaning and the replacement of parts where wear is normal, are not covered by the terms of this guarantee. Visiomed® and its representatives and agents will not in any case be held responsible for any damage and consecutive damages due to the mishandling of this

product. This guarantee is the only valid one at Visiomed[®], any other guarantee (commercial guarantee) except this one will not be taken into account.

IMPORTANT: During the guarantee period if you are dissatisfied with the repairs of this product, please contact the Visiomed[®] customer service.

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EN: *Disposal of Waste Equipment by Users in Private Households in the European Union.*

The symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

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© 06/2016 VM-ZX1
EN: EasyScan™ is a brand registered by Visiomed Group SA. The constructor keeps the right to modify the photographs and characteristics without previous notice. Registered design. The logos and trademarks are the property of their respective owners.