KINECARE[®]

VISIOMED GROUP SA
112, Avenue Kléber
75116 PARIS - France

Made in China

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INSTRUCTIONS FOR USE

CE
0197

This Kinecare® TENS/EMS electrostimulator is a class IIa device that meets all applicable European standards.



TABLE OF CONTENTS

1. PACKAGE CONTENTS	17
2. SAFETY INFORMATION	17
2.1. WARNINGS.....	17
2.2. PRECAUTIONS	18
2.3. GENERAL PRECAUTIONS.....	19
4. INDICATIONS.....	20
5. INSTRUCTIONS BEFORE USE	20
5.1. INSERTING THE BATTERIES	20
5.2. USE OF THE PRODUCT AND TREATMENTS	21
6. CLEANING AND MAINTENANCE.....	23
6.1. CLEANING THE ELECTROSTIMULATOR.....	23
6.2. CLEANING THE PATCHES	23
6.3. ELECTROSTIMULATOR STORAGE CONDITIONS	23
6.4. STORING THE ELECTRODE PATCHES.....	23
6.5. WASTE DISPOSAL	23
7. TROUBLESHOOTING GUIDE	24
8. TECHNICAL FEATURES	25

Dear Customer,
You have just purchased a TENS/EMS electrostimulator Kinecare® from VISIOMED®. In order to use it optimally and efficiently, we recommend that you read these operating instructions carefully.

1. PACKAGE CONTENTS

- 1 electrostimulator
- 1 multi-area self-adhesive patch with protective film (reference VM-MP1)
- 2 batteries

Patches sold separately

VM-AP1
Abdominal patch
170 X 145 mm
6.6 x 5.7 in

VM-MP1
Multizone Patch
237 X 70 mm
9.3 x 2.7 in

VM-SNP1
Patch
Shoulders & Cervical
280 x 90 mm
11.4 x 3.5 in

2. SAFETY INFORMATION

It is important to read all the warnings and precautions in this manual. They are for your safety, to prevent injuries and avoid situations that could damage the device

2.1. WARNINGS

Consult your doctor if you are in any doubt, if you have a serious illness or if you want to know more about using this medical device.

DO NOT USE THIS EQUIPMENT WITH THE FOLLOWING MEDICAL DEVICES:

Using this electrostimulator with other electronic medical devices may cause these devices to malfunction.

Do not use this equipment with the following medical devices:

- An implanted electronic medical device (such as a pacemaker).
- An electronic medical device for maintaining life support (such as a breathing device).
- An implant or metal device (such as a drug delivery system).
- High-frequency surgical equipment.

DO NOT USE THE DEVICE UNDER THE FOLLOWING CONDITIONS:

For hygiene reasons, patches must not be shared with other people. This device is for your own personal use. There is a risk of infection.

- If you are pregnant, ask your doctor or midwife for advice. Do not use on the uterus or abdomen. The effects of TENS on foetal development are not yet known.

THIS DEVICE MUST NOT BE USED IN THE FOLLOWING SITUATIONS:

- Epilepsy (diagnosed or suspected).
- Heart disease (diagnosed or suspected).
- A predisposition to internal bleeding due to impact or injuries.
- In cases of venous or arterial thrombosis, or thrombophlebitis (major arterial circulation problems in the legs).
- Tendonitis or muscular contraction, without consulting your doctor.
- If you have had recent surgery.
- After acute trauma or a fracture. Stimulation can disrupt the healing process.
- A cognitive impairment
- Diabetes (muscle contractions stimulate glucose absorption).

Follow your doctor's instructions.

If you are taking pain medication, ask a doctor before using this device.

THE DEVICE MUST NOT BE USED ON THE FOLLOWING PEOPLE:

- Children or infants. The device has not been evaluated for paediatric use.
- This device is not intended for use by persons (including children) with reduced physical, sensory or mental abilities, or by people who are inexperienced or who have insufficient knowledge, unless in the presence of a person responsible for their safety, under supervision or following the device's operating instructions.

DO NOT USE THIS DEVICE DURING THE FOLLOWING ACTIVITIES:

- In the bath or shower. Water penetration can cause the device to malfunction.
- While sleeping.
- When driving, using machines or similar activities requiring precise control; electrical stimulation can cause short involuntary muscle movements.

DO NOT POSITION THE ELECTRODES ON THE FOLLOWING AREAS:

- Do not apply to the front of the neck (left and right sides, carotid artery), face or head. Stimulation can cause serious adverse reactions affecting your blood pressure or heart rate, and can also cause muscle spasms that can block your airways.
- Do not position the electrode on your chest (heart). Applying an electrical current to this area can cause heart rhythm disorders, which can be fatal.
- Do not apply to open wounds, an inflammation, burns, or skin that is irritated or prone to eczema. Stimulation should not be applied to areas where thrombosis or thrombophlebitis is present as this could facilitate the circulation and lead to an embolism.
- Do not apply the adhesive parts of the patch to your spine. Place them on either side as shown in the image (Chapter 5.2.).
- Do not place the electrodes on your genitals. They may stimulate inappropriate muscles or organs.
- Do not place the electrodes on cancerous lesions.
- Do not position the electrodes on insensitive areas of the body. Burns may occur if the user does not realize that the current is too strong.
- Avoid areas with extensive hair growth (as the patches will not adhere properly).

2.2. PRECAUTIONS

- Do not use for more than 30 minutes per session without seeking medical advice. The long-term effects of electrical stimulation are currently unknown.
- Apply the electrodes to skin that is healthy, clean and dry.
- If skin irritation occurs, discontinue treatment. If the problems continue, see a doctor. Hypersensitivity may occur in some cases.
- Remove the electrodes after use.
- Do not apply ointment or solvent on the skin or electrodes, as this may damage the electrodes and cause a malfunction.
- If you experience side effects, stop using this device and see a doctor.
- See a doctor if your pain has not been diagnosed.
- The Kinecare® electrostimulator treats the symptoms and, as such, suppresses the pain that acts as a protective mechanism. This device is not effective in treating the cause of your pain, including headaches.
- The Kinecare® electrostimulator is not a substitute for pain medications or other pain treatments.
- If your pain does not improve, if it becomes chronic or severe, or continues for more than five days, stop using the device and see your doctor.

- TENS¹ devices do not heal diseases or injuries.

¹TENS: Transcutaneous Electrical Nerve Stimulation. Neuro Stimulation Electrique Transcutanée

Side effects

There are few side effects. These are usually skin irritations (a rash) related to the electrodes. They generally disappear within two hours of the electrostimulation session. See your doctor if they continue for more than 24 hours.

Users may also experience myoclonia (involuntary muscle contractions) when the stimulation is too intense, or sensory disturbances.

If you experience side effects during use such as nausea, fainting, electrode burns, pain or headaches, stop using the device and see your doctor.

You may experience temporary muscle pain after an EMS² session.

²EMS: Electrical Muscle Stimulation. Stimulation Electrique Musculaire

Important information

- Do not bring the electrodes in contact with metal objects during use.
- Do not bend the electrodes.
- Always switch off the device before removing or moving the electrodes.
- Change the electrodes if they are damaged, less adhesive or if you start feeling discomfort during stimulation (uncomfortable pinching or tingling). Do not use a plaster or adhesive tape to secure the electrodes to the skin.
- Do not use electrodes that have expired.
- The patches are intended to be used 15-20 times.
- To store the electrodes, follow the instructions in this manual.
- Ensure that the components are properly connected and that the electrodes are correctly positioned.
- Do not hold the device during use.

2.3. GENERAL PRECAUTIONS

- Read this manual carefully and always follow the treatment instructions.
- Use this device only as intended and as described in this manual.
- The patient is intended operator.
- This device may be used by individuals at home.
- This device is designed for external use.
- This product should not encourage you to self-medicate or adapt your treatment.
- Use and storage: Refer to the technical features (range of room temperature, relative humidity, atmospheric pressure).
- Keep this device somewhere clean and dry.
- Do not expose this device to direct sunlight, water or other liquids, chemical solvents, dust or high temperatures.
- Do not expose this device to electric shocks.
- Protect the electrostimulator from possible impacts.
- Follow the maintenance instructions in this manual.
- Do not service and maintain the device while the device is in use.
- Do not attempt to open the device. If you encounter any problems, contact your dealer.
- Keep out of the reach of children. The device contains small parts that can be swallowed. If swallowed, seek medical advice immediately.
- Check the device before use and stop using it if there is an anomaly, malfunction or in case of discomfort. Prolonged use of a damaged product can cause burns, poor results or a more serious hazard.
- If the device is functioning abnormally, stop using and contact the After-Sales service. Refer to the "Troubleshooting guide" section.
- Do not modify the device or its accessories. This could cause a malfunction and affect the warranty.
- Use this device only with the accessories recommended by the manufacturer.
- Do not modify this equipment without the manufacturer's consent.
- Do not use near electronic monitoring devices (such as heart monitors or electrocardiogram alarms). They may not function properly when used at the same time.
- Do not operate the device within 1 meter of shortwave or microwave medical equipment. Proximity of such equipment to this device may deteriorate the performance of the device.
- Portable RF communication devices (including peripheral devices such as antenna cables and external antennas) must not be used within 30 cm of the device. They could deteriorate the performance of the device.
- The use of accessories or devices other than those mentioned or supplied by this equipment's

manufacturer could result in increased electromagnetic emissions or a decrease in the electromagnetic immunity of this equipment and increase the risk of a malfunction.

3. INFORMATION - TRANSCUTANEOUS ELECTRICAL NERVE STIMULATION

The Kinecare® TENS/EMS electrostimulator is a device to relieve adults of acute post-surgical pain or chronic mild to moderate pain. It works on the principle of transcutaneous electrical nerve stimulation (TENS), which helps relieve pain and soothe muscle tension. This device also helps promote strengthen muscle mass. The device has been specially designed for use in the home.

Transcutaneous Electrical Nerve Stimulation (TENS) is a non-medical pain-relief treatment with no side effects used by doctors and physiotherapists for more than 30 years. Electrodes are placed on the relevant areas of the body and send a mild electric current to the nerves. This stimulation causes the body to produce and diffuse natural painkillers (endorphin and enkephalin).

In case of pain or muscle tension, mini electrical pulses in the nerve tissues block pain signals sent to the nervous system and trigger the release of endorphins.

To strengthen your muscle mass, just select the electrical muscle stimulation program (EMS). The Kinecare® TENS/EMS electrostimulator's electrodes, placed on the area to be stimulated, send electrical pulses causing the muscles to contract, thus strengthening them.

The device is operated directly on the Kinecare® TENS/EMS electrostimulator. A diagram of the human body with treatment areas helps you position the electrodes correctly (see section 5.2.).

4. INDICATIONS


- Relief from pain and muscle tension (TENS: 3 pre-set programs).
- Strengthening muscle mass (EMS: 1 pre-set program).
- Relaxation and relief from muscle tension (MASSAGE program).
- Intended for temporary pain relief.
- Designed to stimulate healthy muscles.

Intended use:

The device serves the transcutaneous electrical stimulation of nerves and muscles in humans (muscle and joint pain)

5. INSTRUCTIONS BEFORE USE

5.1. INSERTING THE BATTERIES

1. Remove the cover of the battery compartment using a sharp instrument, such as the tip of a pen, or a flat screwdriver.

2. Insert the batteries AAA into the compartment, following the + and - markings

3. Reposition the battery compartment cover.

If the stimulation intensity is too low, if the power of the device is low, or if the device does not respond, replace the batteries.

Caution:

- Remove the batteries if the device is not used for long periods.
- Keep batteries out of the reach of children.
- Dispose of batteries safely in accordance with local regulations.
- Do not mix old and new batteries or different types of batteries.
- If the batteries leak and come into contact with your skin or eyes, rinse immediately and abundantly and see a doctor.

NB: The approximate lifespan of the batteries is 30 sessions.

5.2. USE OF THE PRODUCT AND TREATMENTS

Before positioning the patch, check that the skin is free of wounds or rashes and that it is healthy.

- Transcutaneous Electrical Nerve Stimulation (TENS) programs: Do not position the patch on skin that is irritated or that has open wounds, a rash or cancerous lesions.

- Electrical Muscle Stimulation (EMS) programs: Position the patch only on muscles that are healthy and intact.

TENS program

A usual transcutaneous nerve stimulation program uses high-frequency stimulation. The purpose of this stimulation is to block the pain signal to the nervous system. Typically, people use conventional transcutaneous nerve stimulation for a session lasting around 30 minutes. It is possible to use these programs for different durations and at different frequencies as needed. It is normal to feel the effects of this program more intensely during the stimulation and for the effects to fade after the session. To avoid dependence, it is advisable to start the program on a low intensity and to increase intensity during the session, while avoiding pain or discomfort.

EMS program

These programs, commonly known as "neuromuscular electrical stimulation" (NMES) or "electromyostimulation", produce pulses causing the muscle to contract. They are used in particular to strengthen muscle mass. Sessions should be carried out at rest (sitting or lying).

When used for the first time, electrical stimulation can cause an unusual and strange sensation. To avoid dependence it is advisable to start the program on a low intensity and to increase intensity during the session. The stimulation should remain comfortable and the muscle contractions should not be tiring or painful. If the program is used after an intense physical exercise or after exertion, it is advisable to choose a low intensity to limit muscle fatigue.

Step 1: Cleaning the skin

Clean your skin using a damp cloth then dry. Any oil, cream or lotion, dust or other material on the skin may prevent the patch from adhering. Avoid areas with extensive hair growth.

Step 2: Clip the patch onto the electrostimulator

Before use, check the expiry date and condition of the electrodes (deterioration, wear).

Before switching on the device, attach the patch to the device using the press studs.

Step 3: Positioning the patch

Remove the protective film from the electrodes. Keep this film to reposition on the patch after the session.

Place the electrodes on the treatment area. See the patch positioning diagrams detailed below in this manual. Do not move the electrodes while the product is in use. Stop the device before moving the patch.

Caution

As stated in the warnings "DO NOT POSITION THE ELECTRODES ON THE FOLLOWING AREAS":

The electrodes must never be placed on any of the following areas: the head, the carotid artery, on either side of the neck, at the front of the neck or on the chest (heart area). Do not place the electrode patch on skin that is inflamed or that has open wounds, a rash or cancerous lesions or on irritated skin.

The effectiveness of your sessions will depend in part on how you place the patch on the skin. It is recommended to place it on or around the area to be treated.

Before positioning the patches, ensure that the area is sufficiently sensitive. If you are unsure, do not hesitate to move the patch for greater sensitivity. (Stop the stimulation by pressing on the button  before moving the device.)

In case of proven effectiveness in a specific area, remember the location of the electrodes to place them there again during a later session.

If you experience discomfort (unpleasant contractions, particularly during the electrical nerve stimulation programs in bursts), it is advisable to move the patch. If, during the electrical muscle stimulation programs (intended in particular to contract the muscle) no contraction occurs, either the intensity is too low or the patch is not positioned properly.

Program	Mode	Pad Application
1	TENS	
2	TENS	
3	TENS	
4	EMS	
5	Massage	

NB: The neck & shoulder patch and the abdomen patch are sold separately.

The Massage program can be used on the back areas recommended for other programs.

Step 4: Switching the electrostimulator on

Press on the button  to switch on the device.

Step 5: Select the treatment program

The wireless TENS / EMS electrostimulator has five treatment programs. Press on the button  to select the treatment program you need.

The LED showing the selected program from 1 to 5 lights up.

Step 6: Start the program and adjust the intensity.

Press on the + button to start a program. The LED indicator stays on and shows the program number.

Once the program is selected, the intensity can be changed using the + and - buttons. The maximum intensity level is 15.

NB:

You can press on the button  to stop the program. The electrostimulator instantly sets the intensity level to zero.

CAUTION:

- If the electrodes are not positioned properly on the skin or if the device is not connected to the electrodes and the intensity level is 1 or greater than 1, the device output will automatically shut off.
- If the stimulation levels are or become uncomfortable, reduce the intensity to a comfortable level. Contact your doctor if the problems continue.
- If the pain does not improve and you are suffering after excessive use, do not treat these areas for two days. If this problem continues, reduce the treatment duration and intensity.
- If you experience pain, dizziness, discomfort or nausea, contact your doctor.
- Treatment at the maximum intensity should be used with caution. Do not exceed your comfort level (or threshold).
- Do not switch on the device if the electrodes are not positioned on your body.
- Never remove the patch from your skin during the session. Always switch off the power before removing the patch.

NB:

- It is possible to carry out a session with your clothes on. Simply place the electrodes under your clothes.
- You may need help from another person to place the electrodes in certain areas.

Step 7: Switching the electrostimulator off

There are several possibilities:

1. Press the button  and hold for two seconds to switch the power off.
2. The device switches off automatically after 3 minutes of inactivity.

6. CLEANING AND MAINTENANCE

6.1. CLEANING THE ELECTROSTIMULATOR

1. Switch off the device.
2. Clean the device after use with a soft, slightly damp cloth. Wipe gently. Do not use chemicals (e.g. benzene). Do not immerse the device in any type of liquid. Never allow water to enter the battery compartment or come into contact with the snap fasteners.

NB:

This device and its accessories do not require sterilization.

6.2. CLEANING THE PATCHES

1. Switch off the device
2. Remove the patches from the device.
3. The patches should be cleaned when the adhesive surface becomes dirty and/or does not adhere properly to your skin. Gently wash the patch with your fingers under cold water for several seconds. Do not use a sponge, cloth or sharp object. Do not use a detergent, chemical or soap.
4. Dry the outer side of the patches with a soft cloth and leave the adhesive surface to air dry completely.

NB:

- Clean the patches after each use.
- The lifespan of the patches may vary depending on how often they are washed, the condition of your skin and how they are stored.
- The patches are intended to be used 15-20 times.
- If the patch no longer adheres to your skin or is damaged, replace it with a new patch.
- Only use KINECARE® patches, which are specially designed for this purpose.

6.3. ELECTROSTIMULATOR STORAGE CONDITIONS

Place the electrostimulator, the electrodes and the manual in their original packaging. Keep it in a clean place and away from moisture. See storage conditions in the technical specifications

NB:

Remove the batteries if the device is not used for long periods.

6.4. STORING THE ELECTRODE PATCHES

Place the clean, dry patches on the plastic film and store them in the original packaging between each use.

6.5. WASTE DISPOSAL

 Dispose of the device in accordance with Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE). Contact your local authorities for waste disposal. Patch: This product is not hazardous waste. Please dispose of the product in accordance with local regulations.

7. TROUBLESHOOTING GUIDE

Problem	Possible cause	Solution
The device does not switch on	The batteries are discharged	Put new batteries in the device
	The batteries are not properly installed	Install the batteries properly, following the correct polarity
Treatment appears to be different or less pleasant than previously	The patch is positioned incorrectly (wrong place, wrongly applied, etc.)	Switch off the device and reposition the patch on your skin
	The intensity is too low or too high	Change the intensity using the + or - buttons
	The electrodes' expiry date is exceeded	Replace the patch
	The patch is deteriorated or worn	Replace the patch
	The patch is dirty	Clean the patch, if the problem continues, replace the patch
Unpleasant sensation when touching the electrode or device fasteners	The session was not stopped correctly (the on/off button was not pressed before removing the device).	Always switch off the device at the end of the session before removing the electrodes. Press the STOP button for 2 seconds.
Your skin becomes red and/or you feel pain	The patches are always used in the same place	Position the patches differently
	The patch does not adhere properly to your skin	Make sure the patches are positioned correctly on your skin
	The patch is dirty	Clean the patch. If the problem continues, replace the patch
	The surface of the patch is damaged	Replace the patch

8. TECHNICAL FEATURES

Brand	KINECARE®
Model	VM-WT1 Tens/EMS (wireless)
Units	1 device and 1 patch (reference VM-MP1)
Programs	5 programs (3 TENS, 1 EMS, 1 Massage)
Classification	Medical device Class IIa
Intensity levels	0 to 15
Power supply	2 AAA batteries
Output voltage	45V (charged) at 500 ohm
Operating conditions	5°C - 40°C / 41°F - 104°F with a relative humidity of 30% to 75% and atmospheric pressure of 700 to 1060 hPa.
Storage and transport conditions	-10°C - 55°C / 14°F - 131°F with a relative humidity of 10% to 90% and atmospheric pressure of 700 to 1060 hPa.
Electrostimulator dimensions	71mm(L)x55mm(w)x18.5mm(H) 2.79in(L)x2.16in(l)x0.72(H)
Weight	38g / 0,08lb (without batteries)
Lifespan of the batteries	Approx. 30 uses (if used 30 min/day at half the maximum intensity with new batteries on program 1)
Service life of device	2 years
IP rating	IP22
Waveform	Biphasic square wave
Pulse width	30 µs-250 µs
Frequency	2Hz-290 Hz
Treatment time	30 minutes
Electrode dimensions	210mm(L)x68mm(w) 8.26in(L)x2.67in(l)

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

VM-AP1
Buikpatch
170 X 145 mm
6.6 x 5.7 in

VM-MP1
Multizone Patch
237 X 70 mm
9.3 x 2.7 in

VM-SNP1
lap
Schouders & Cervicale
280 x 90 mm
11.4 x 3.5 in

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ELECTRO-
MAGNETIC
COMPATIBILITY
DATA

EMC data

Important information regarding electromagnetic compatibility (emc).

Table 1

Declaration - electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	Complies
Harmonic emissions IEC 61000-3-2	Not applicable	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Table 2

Declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Compliance	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315°	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
	0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles		
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

Table 3

Declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3V 0.15 MHz to 80MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of device, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ 150 KHz to 80 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which device is used exceeds the applicable RF compliance level above, device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating device.
Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4

Recommended separation distances between portable and mobile RF communications equipment and device			
Device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and device, as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	0.15 MHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	80 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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EN: Visiomed* will repair or replace this product free of charge in the case of defective parts or manufacturing defects, in accordance with the conditions mentioned below as follows:

DURATION: 24 MONTHS RETURN TO WORKSHOP

LIMITS AND EXCLUSIONS: This warranty applies only to the original final purchaser. A purchase invoice, or other proof of purchase, and this warranty card will be requested for after-sales service in accordance with this warranty. This warranty card will not be extended to anyone other than the original final purchaser. This warranty is voided if the serial numbers on the product are changed, replaced, illegible or missing, or if repairs have been attempted by any unauthorized service, including the user. This warranty covers only material or parts defects that occur during normal use of the product. It does not cover damages caused during shipment or transportation of the device, caused by repairs carried out by a distributor, by alterations made, by connecting equipment that has not been approved by Visiomed, or by any use failing to comply with the operating instructions. In addition, this warranty does not cover damage from falls, improper handling, improper installation, damage from fire, flooding, lightning, or any other natural disaster. This warranty does not cover the packaging of the equipment, accessories, defects in appearance due to the commercial display of the product, showroom, sales area, demonstration, etc. Normal maintenance, cleaning and the replacement of parts undergoing normal wear are not covered by this warranty.

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EN: *Disposal of Waste Equipment by Users in Private Households in the European Union.*

For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

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EN: *Kinecare® is a brand registered by Visiomed Group SA. The constructor keeps the right to modify the photographs and characteristics without previous notice. Registered design. The logos and trademarks are the property of their respective owners.*

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© 10/2017 KINECARE